



## The Classic

### Professional Design

LED display, easy operation, stainless steel water tank, high density bacteria filter.

### Leading Technology

Automatic sterilization, automatic drying, keep-temperature function, adjustable sterilization time, accurate temperature control, safety protection design, monitoring system.

Program	121°C Solid	121°C Universal	134°C Solid	134°C Universal	Cotton	BD/Heix Test	Vacuum Test
Vacuum times	1	3	1	3	3	3	3
Temperature (°C)	121	121	134	134	134	134	/
Pressure (KPA)	110	110	210	210	210	210	-80
Time (Min.)	20	20	4	4	4	3.5	15

Class B



## The Competitive

### Professional Design

Filtered water can also be used, automatic water filling, automatic water-release, LED display, easy operation, stainless steel water tanks, high density bacteria filter.

### Leading Technology

Steam generated inside chamber, automatic sterilization & dry, keep-temperature function adjustable sterilization time, accurate temperature control, safety protection design, easy maintenance.

Program	121°C Solid	121°C Universal	134°C Solid	134°C Universal	Cotton	BD/Heix Test	Vacuum Test
Vacuum times	1	3	1	3	3	3	3
Temperature (°C)	121	121	134	134	134	134	/
Pressure (KPA)	110	110	210	210	210	210	-80
Time (Min.)	20	20	4	4	4	3.5	15

Class B



### YS-8L

Capacity	8L
Chamber Size	φ17*32cm
Display	LED
Lock	Electrical/Mechanical
Product Size	53.5*36.5*32cm
G.W./N.W.	40/34 KGS
Power	220V/50HZ 1.7KW
Output	RS232
Printer	thermal/ stylus

### YS-12L

Capacity	12L
Chamber Size	φ20*36cm
Display	LED
Lock	Electrical/Mechanical
Product Size	66*44.5*39cm
G.W./N.W.	53/45 KGS
Power	220V/50HZ 1.7KW
Output	RS232
Printer	thermal/ stylus

### YS-18L

Capacity	18L
Chamber Size	φ25*35cm
Display	LED
Lock	Electrical/Mechanical
Product Size	66*44.5*39cm
G.W./N.W.	56/48 KGS
Power	220V/50HZ 1.7KW
Output	RS232
Printer	thermal/ stylus

### YS-22L

Capacity	22L
Chamber Size	φ25*45cm
Display	LED
Lock	Electrical/Mechanical
Product Size	70*44.5*39cm
G.W./N.W.	58/50 KGS
Power	220V/50HZ 1.7KW
Output	RS232
Printer	thermal/ stylus

### YS-8L-II

Capacity	8L
Chamber Size	φ17*32cm
Display	LED
Lock	Electrical/Mechanical
Product Size	53.5*36.5*32cm
G.W./N.W.	34/28 KGS
Power	220V/50HZ 1.7KW
Output	RS232
Printer	thermal/ stylus

### YS-12L-II

Capacity	12L
Chamber Size	φ20*36cm
Display	LED
Lock	Electrical/Mechanical
Product Size	66*44.5*39cm
G.W./N.W.	46/37 KGS
Power	220V/50HZ 1.7KW
Output	RS232
Printer	thermal/ stylus

### YS-18L-II

Capacity	18L
Chamber Size	φ25*35cm
Display	LED
Lock	Electrical/Mechanical
Product Size	66*44.5*39cm
G.W./N.W.	49/40 KGS
Power	220V/50HZ 1.7KW
Output	RS232
Printer	thermal/ stylus

### YS-22L-II

Capacity	22L
Chamber Size	φ25*45cm
Display	LED
Lock	Electrical/Mechanical
Product Size	70*44.5*39cm
G.W./N.W.	52/43 KGS
Power	220V/50HZ 1.7KW
Output	RS232
Printer	thermal/ stylus



yeson@yeson-medicine.com  
www.yeson-medicine.com

yeson@yeson-medicine.com  
www.yeson-medicine.com

