

Electrosurgical Unit  
**ZEUS-100/80**



**Output Specification**

Mode	ZEUS-100	ZEUS-80
Pure Cut	100W /300Ω	80W /300Ω
Blend Cut	100W /300Ω	80W /300Ω
Contact Coagulation	60W /300Ω	50W /300Ω
Bipolar Coagulation	80W /100Ω	70W /100Ω



- Main Frequency : 1.6MHz
- Power Source : AC120V or AC230V, 50Hz / 60Hz
- Protection Class : Class 1, Type BF
- Dimension(WxDxHmm) : 270×310×110
- Weight : 6.0Kg

**Standard Accessories**

Code No.	Name of Item
FS01-01R	Single Foot Switch
HD02-01D	Disposable Twin Button Handle
PL02-05R	Silicon Patient Plate(165X90mm)
EL01-02D	Electrode(Knife 2.4*70mm)
EL02-02D	Electrode(Ball 5mm)
EL04-02D	Electrode(Needle 2.4*70mm)
EL04-04D	Electrode(Needle Angled 2.4*70mm)



FS01-01R



HD02-01D



PL02-05R




EL01-02D / EL02-02D  
EL04-02D / EL04-04D

Electrosurgical Unit  
**ZERO 50**



**Output Specification**

Mode	ZERO50
Cut	50W /300Ω
Coagulation	40W /300Ω



- Main Frequency : 1.6MHz
- Power Source : AC120V or AC230V, 50Hz / 60Hz
- Protection Class : Class 1, Type BF
- Dimension(WxDxHmm) : 200×220×90
- Weight : 2.6Kg

**Standard Accessories**

Code No.	Name of Item
HD02-01D	Disposable Twin Button Handle
PL02-03R	Silicon Patient Plate(110X80mm)



HD02-01D



PL02-03R