




RP11A

FLEXIBLE VIDEO RHINOLARYNGOSCOPE



160g

Ultra-light handle, more stable grip, reduce the fatigue of the doctor, and make the examination smoother.



1280x720P

3.2mm high-definition video endoscope, ultra-thin insertion, full of details

A variety of connection methods are suitable for various diagnosis and treatment environments.



Computer



Pad



Processor

Optical system

- Field of view 100°
- Depth of field 3-150mm
- Direction of view 0°

Insertion portion

- Insertion tube OD 3.2mm/3.8mm
- Angulation range U 130° / D 130°
- Working length 330mm±10%
- Extension cable USB、REMO
- Max. resolution 1280x720P



weseemed
medical technology

ZHONGSHAN WESEE MEDITECH CO., LTD.

6F, Block C, 2Jianshuo Str. HoChung Village,
Nanlang, Zhongshan 528451 P. R. China

TEL:(+86)0760-86960005

Email:info@weseemed.com

www.weseemed.com

